

Product datasheet for **GA203214**

Mouse **Pik3c2g** activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pik3c2g CRISPRa kit - CRISPR gene activation of mouse phosphatidylinositol-4-phosphate 3-kinase catalytic subunit type 2 gamma
Symbol:	Pik3c2g
UniProt ID:	<u>Q70167</u>
Synonyms:	C80387; PI3K-C2-gamma
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)