

Product datasheet for **GA203215**

Mouse *Pik3ca* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Pik3ca</i> CRISPRa kit - CRISPR gene activation of mouse phosphatidylinositol-4,5-bisphosphate 3-kinase catalytic subunit alpha
Symbol:	<i>Pik3ca</i>
UniProt ID:	P42337
Synonyms:	6330412C24Rik; caPI3K; p110; p110alpha
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)