

## Product datasheet for **GA203217**

### Mouse **Pik3r1** activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Pik3r1 CRISPRa kit - CRISPR gene activation of mouse phosphoinositide-3-kinase regulatory subunit 1
Symbol:	Pik3r1
UniProt ID:	<a href="#">P26450</a>
Synonyms:	p50alpha; p55alpha; p85alpha; PI3K
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)