

Product datasheet for **GA203219**

Mouse **Pik3r3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Pik3r3 CRISPRa kit - CRISPR gene activation of mouse phosphoinositide-3-kinase regulatory subunit 3
Symbol:	Pik3r3
UniProt ID:	Q64143
Synonyms:	AA414954; p55pik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)