

Product datasheet for **GA203220**

Mouse **Pikfyve** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Pikfyve CRISPRa kit - CRISPR gene activation of mouse phosphoinositide kinase, FYVE type zinc finger containing
Symbol:	Pikfyve
UniProt ID:	<u>Q9Z1T6</u>
Synonyms:	5230400C17Rik; p235; Pip5k; Pip5k3; Pipk5k3; PipkIII
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)