

Product datasheet for **GA203224**

Mouse Pip5k1c activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pip5k1c CRISPRa kit - CRISPR gene activation of mouse phosphatidylinositol-4-phosphate 5-kinase, type 1 gamma
Symbol:	Pip5k1c
UniProt ID:	Q70161
Synonyms:	A1115456; A1835305; Pip5klgamma
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)