

## Product datasheet for **GA203225**

### Mouse Pip4k2a activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Pip4k2a CRISPRa kit - CRISPR gene activation of mouse phosphatidylinositol-5-phosphate 4-kinase, type II, alpha
Symbol:	Pip4k2a
UniProt ID:	<u><a href="#">Q70172</a></u>
Synonyms:	AW742916; Pip5k2a
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)