

Product datasheet for **GA203230**

Mouse Pirb activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pirb CRISPRa kit - CRISPR gene activation of mouse paired Ig-like receptor B
Symbol:	Pirb
Components:	GA203230G1, Pirb gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA203230G2, Pirb gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA203230G3, Pirb gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P97484</u>
Synonyms:	Gp91; Lilrb3; LIR-3; PIR-B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)