

## Product datasheet for **GA203236**

### Mouse Pitx3 activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | Pitx3 CRISPRa kit - CRISPR gene activation of mouse paired-like homeodomain transcription factor 3 |
| Symbol:       | Pitx3  |
| UniProt ID:   | <a href="#">O35160</a>   |
| Synonyms:     | ak; Ptx3   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                            |

#### Product images:



[View online »](#)