

Product datasheet for **GA203248**

Mouse Prkd1 activation kit by CRISPRa

Product data:

| | |
|---------------|-------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits |
| Description: | Prkd1 CRISPRa kit - CRISPR gene activation of mouse protein kinase D1 |
| Symbol: | Prkd1 |
| UniProt ID: | Q62101 |
| Synonyms: | Pkcm; PKD; PKD1; Prkcm |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)