

Product datasheet for **GA203251**

Mouse Pkd2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pkd2 CRISPRa kit - CRISPR gene activation of mouse polycystin 2, transient receptor potential cation channel
Symbol:	Pkd2
UniProt ID:	Q35245
Synonyms:	C030034P18Rik; PC2; TRPP2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)