

Product datasheet for **GA203257**

Mouse Pknox1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pknox1 CRISPRa kit - CRISPR gene activation of mouse Pbx/knotted 1 homeobox
Symbol:	Pknox1
Components:	GA203257G1, Pknox1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA203257G2, Pknox1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA203257G3, Pknox1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O70477</u>
Synonyms:	D17Wsu76e; PREP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)