

Product datasheet for **GA203281**

Mouse *Plcd4* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Plcd4 CRISPRa kit - CRISPR gene activation of mouse phospholipase C, delta 4
Symbol:	Plcd4
UniProt ID:	Q8K3R3
Synonyms:	4921507K24Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)