

Product datasheet for **GA203287**

Mouse Prl2c3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Prl2c3 CRISPRa kit – CRISPR gene activation of mouse prolactin family 2, subfamily c, member 3
Symbol:	Prl2c3
Components:	GA203287G1, Prl2c3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA203287G2, Prl2c3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA203287G3, Prl2c3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P04768</u>
Synonyms:	Ghd1; MRP-2; MRP-3; mrp/plf3; Mrpplf3; PLF-2; PLF-3; Plf2; Prl2c4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

