

Product datasheet for **GA203294**

Mouse Plp2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Plp2 CRISPRa kit - CRISPR gene activation of mouse proteolipid protein 2
Symbol:	Plp2
Components:	GA203294G1, Plp2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA203294G2, Plp2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA203294G3, Plp2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9R1Q7</u>
Synonyms:	mIMA4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)