

Product datasheet for **GA203301**

Mouse *Plxna3* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Plxna3 CRISPRa kit - CRISPR gene activation of mouse plexin A3
Symbol:	Plxna3
UniProt ID:	P70208
Synonyms:	mKIAA4078; PlexA3; Plxa3; Plxn3; Plxn4; SEX
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)