

Product datasheet for **GA203344**

Mouse Ppef2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ppef2 CRISPRa kit - CRISPR gene activation of mouse protein phosphatase, EF hand calcium-binding domain 2
Symbol:	Ppef2
Components:	GA203344G1, Ppef2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA203344G2, Ppef2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA203344G3, Ppef2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q35385</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)