

## Product datasheet for **GA203349**

### Mouse Lgals3bp activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Lgals3bp CRISPRa kit - CRISPR gene activation of mouse lectin, galactoside-binding, soluble, 3 binding protein
Symbol:	Lgals3bp
Components:	GA203349G1, Lgals3bp gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA203349G2, Lgals3bp gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA203349G3, Lgals3bp gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q07797</a></u>
Synonyms:	90K; CyCAP; MAC-2BP; Ppicap; Tango10b
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)