

Product datasheet for **GA203362**

Mouse Ppp3cb activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Ppp3cb CRISPRa kit - CRISPR gene activation of mouse protein phosphatase 3, catalytic subunit, beta isoform |
| Symbol: | Ppp3cb |
| UniProt ID: | P48453 |
| Synonyms: | 1110063J16Rik; Calnb; Cnab; CnAbeta |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)