

Product datasheet for **GA203364**

Mouse Ppp3r1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ppp3r1 CRISPRa kit - CRISPR gene activation of mouse protein phosphatase 3, regulatory subunit B, alpha isoform (calcineurin B, type I)
Symbol:	Ppp3r1
UniProt ID:	Q63810
Synonyms:	CaNB1; Cnb1; MCIP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)