

Product datasheet for **GA203371**

Mouse Mob4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mob4 CRISPRa kit - CRISPR gene activation of mouse MOB family member 4, phocein
Symbol:	Mob4
UniProt ID:	<u>Q6PEB6</u>
Synonyms:	2610109B12Rik; Mobkl3; Phocn; Prei3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)