

Product datasheet for **GA203389**

Mouse Prl4a1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Prl4a1 CRISPRa kit - CRISPR gene activation of mouse prolactin family 4, subfamily a, member 1
Symbol:	Prl4a1
Components:	GA203389G1, Prl4a1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA203389G2, Prl4a1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA203389G3, Prl4a1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q35256</u>
Synonyms:	AA409771; Ghdl23; PLP-A; Prlpa
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)