

Product datasheet for **GA203424**

Mouse Psm3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Psm3 CRISPRa kit - CRISPR gene activation of mouse proteasome (prosome, macropain) subunit, alpha type 3
Symbol:	Psm3
UniProt ID:	O70435
Synonyms:	Lmpc8
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)