

Product datasheet for **GA203430**

Mouse Psmc3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Psmc3 CRISPRa kit - CRISPR gene activation of mouse proteasome (prosome, macropain) 26S subunit, ATPase 3
Symbol:	Psmc3
UniProt ID:	<u>O88685</u>
Synonyms:	TBP-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)