

## Product datasheet for **GA203431**

### Mouse Psmc3ip activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Psmc3ip CRISPRa kit - CRISPR gene activation of mouse proteasome (prosome, macropain) 26S subunit, ATPase 3, interacting protein
Symbol:	Psmc3ip
UniProt ID:	<u><a href="#">Q35047</a></u>
Synonyms:	C79099; GT198; HOP2; Tbpip
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)