

## Product datasheet for **GA203434**

### Mouse Psme3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Psme3 CRISPRa kit - CRISPR gene activation of mouse proteasome (prosome, macropain) activator subunit 3 (PA28 gamma, Ki)
Symbol:	Psme3
UniProt ID:	<a href="#">P61290</a>
Synonyms:	AA410043; AU020960; Ki; pa28g; PA28gamma; REGgamma
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)