

Product datasheet for **GA203451**

Mouse *Ptger1* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Ptger1</i> CRISPRa kit - CRISPR gene activation of mouse prostaglandin E receptor 1 (subtype EP1)
Symbol:	<i>Ptger1</i>
UniProt ID:	P35375
Synonyms:	42kDa; EP1; <i>Ptgerep1</i>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)