

## Product datasheet for **GA203453**

### Mouse *Ptger3* activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	<i>Ptger3</i> CRISPRa kit - CRISPR gene activation of mouse prostaglandin E receptor 3 (subtype EP3)
Symbol:	<i>Ptger3</i>
Synonyms:	EP3; Pgerep3; <i>Ptgerep3</i>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)