

Product datasheet for **GA203454**

Mouse *Ptger4* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Ptger4</i> CRISPRa kit - CRISPR gene activation of mouse prostaglandin E receptor 4 (subtype EP4)
Symbol:	<i>Ptger4</i>
UniProt ID:	P32240
Synonyms:	EP4; <i>Ptgerep4</i>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)