

Product datasheet for **GA203457**

Mouse Ptgir activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ptgir CRISPRa kit - CRISPR gene activation of mouse prostaglandin I receptor (IP)
Symbol:	Ptgir
UniProt ID:	<u>P43252</u>
Synonyms:	IP; PGI2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)