

Product datasheet for **GA203473**

Mouse Ptpn11 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ptpn11 CRISPRa kit - CRISPR gene activation of mouse protein tyrosine phosphatase, non-receptor type 11
Symbol:	Ptpn11
UniProt ID:	P35235
Synonyms:	2700084A17Rik; AW536184; PTP1D; PTP2C; SAP-2; SH-PTP2; SH-PTP3; SHP-2; Shp2; Syp
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)