

Product datasheet for **GA203480**

Mouse Ptpn4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ptpn4 CRISPRa kit - CRISPR gene activation of mouse protein tyrosine phosphatase, non-receptor type 4
Symbol:	Ptpn4
UniProt ID:	Q9WU22
Synonyms:	hPTP-MEG; Ptn4; PTPMEG; TEP; TEP/mPTPMEG
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)