

Product datasheet for **GA203487**

Mouse Ptprcap activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ptprcap CRISPRa kit - CRISPR gene activation of mouse protein tyrosine phosphatase, receptor type, C polypeptide-associated protein
Symbol:	Ptprcap
UniProt ID:	Q64697
Synonyms:	A1605264; CD-45-AP; LPAP; LSM-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)