

Product datasheet for GA203520

Mouse Qk activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Qk CRISPRa kit – CRISPR gene activation of mouse quaking
Symbol:	Qk
Components:	GA203520G1, Qk gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA203520G2, Qk gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA203520G3, Qk gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9QYS9</u>
Synonyms:	I110003F05Rik; I500005P18; I(17)-1W1s; I17W1s1; Qkl
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

