


## Product datasheet for **GA203539**

### Mouse Dennd5a activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | Dennd5a CRISPRa kit - CRISPR gene activation of mouse DENN/MADD domain containing 5A                                                                                                                                                                                                                             |
| Symbol:       | Dennd5a                                                                                                                                                                                                                                                                                                          |
| Components:   | GA203539G1, Dennd5a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA203539G2, Dennd5a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA203539G3, Dennd5a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u>Q6PAL8</u>                                                                                                                                                                                                                                                                                                    |
| Synonyms:     | 1500012B19Rik; AJ245569; AW496491; mKIAA1091; ORF37; Rab6ip1                                                                                                                                                                                                                                                     |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)