

Product datasheet for **GA203560**

Mouse Rai1 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Rai1 CRISPRa kit - CRISPR gene activation of mouse retinoic acid induced 1 |
| Symbol: | Rai1 |
| UniProt ID: | <u>Q61818</u> |
| Synonyms: | Gt1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)