

Product datasheet for **GA203570**

Mouse *Rasa3* activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Rasa3 CRISPRa kit - CRISPR gene activation of mouse RAS p21 protein activator 3 |
| Symbol: | Rasa3 |
| UniProt ID: | Q60790 |
| Synonyms: | A1326412; C86362; E130011G04; GAPIII; hlb381; scat |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)