

Product datasheet for **GA203642**

Mouse Ripk1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ripk1 CRISPRa kit - CRISPR gene activation of mouse receptor (TNFRSF)-interacting serine-threonine kinase 1
Symbol:	Ripk1
UniProt ID:	<u>Q60855</u>
Synonyms:	D330015H01Rik; Rinp; RIP; Rip1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)