

Product datasheet for **GA203646**

Mouse Uri1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Uri1 CRISPRa kit - CRISPR gene activation of mouse URI1, prefoldin-like chaperone
Symbol:	Uri1
UniProt ID:	Q3TLD5
Synonyms:	C80913; NNX3; Rmp
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)