


## Product datasheet for **GA203665**

### Mouse Snord73b activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                     |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                        |
| Description:  | Rnu73b CRISPRa kit - CRISPR gene activation of mouse U73B small nuclear RNA                                                                                                                                                                                                                                         |
| Symbol:       | Snord73b                                                                                                                                                                                                                                                                                                            |
| Components:   | GA203665G1, Snord73b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA203665G2, Snord73b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA203665G3, Snord73b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                             |

#### Product images:



[View online »](#)