

## Product datasheet for **GA203707**

### Mouse Rps14 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                            |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| Description:  | Rps14 CRISPRa kit - CRISPR gene activation of mouse ribosomal protein S14                                                                                                                                                                                                                                  |
| Symbol:       | Rps14                                                                                                                                                                                                                                                                                                      |
| Components:   | GA203707G1, Rps14 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA203707G2, Rps14 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA203707G3, Rps14 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">P62264</a>                                                                                                                                                                                                                                                                                     |
| Synonyms:     | 2600014J02Rik; AL023078                                                                                                                                                                                                                                                                                    |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:



[View online »](#)