

## Product datasheet for **GA203714**

### Mouse Rps3a1 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | Rps3a1 CRISPRa kit - CRISPR gene activation of mouse ribosomal protein S3A1 |
| Symbol:       | Rps3a1  |
| UniProt ID:   | <a href="#">P97351</a>  |
| Synonyms:     | Rps3a   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)     |

#### Product images:



[View online >](#)