

Product datasheet for **GA203747**

Mouse **S100a4** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	S100a4 CRISPRa kit - CRISPR gene activation of mouse S100 calcium binding protein A4
Symbol:	S100a4
UniProt ID:	P07091
Synonyms:	18A2; 42a; Cap1; FSp1; metastasin; Mts1; PeL98; pk9a
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)