

Product datasheet for **GA203750**

Mouse Calprotectin (S100a8) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	S100a8 CRISPRa kit - CRISPR gene activation of mouse S100 calcium binding protein A8 (calgranulin A)
Symbol:	S100a8
UniProt ID:	P27005
Synonyms:	60B8Ag; AI323541; B8Ag; CAg; CFAg; CP-10; MRP8; p8
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)