

Product datasheet for **GA203751**

Mouse Calprotectin (S100a9) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	S100a9 CRISPRa kit - CRISPR gene activation of mouse S100 calcium binding protein A9 (calgranulin B)
Symbol:	S100a9
UniProt ID:	P31725
Synonyms:	60B8Ag; AW546964; BEE22; Cagb; GAGB; L1Ag; MRP14; p14
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)