

Product datasheet for **GA203758**

Mouse Sag activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sag CRISPRa kit - CRISPR gene activation of mouse S-antigen, retina and pineal gland (arrestin)
Symbol:	Sag
UniProt ID:	P20443
Synonyms:	A930001K18Rik; Arr1; arrestin; Irbp
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)