

Product datasheet for **GA203763**

Mouse Sf3a2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sf3a2 CRISPRa kit - CRISPR gene activation of mouse splicing factor 3a, subunit 2
Symbol:	Sf3a2
Components:	GA203763G1, Sf3a2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA203763G2, Sf3a2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA203763G3, Sf3a2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	66kDa; AW047242; PRP11; Sap62; SFA66
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)