

Product datasheet for **GA203779**

Mouse *Scn10a* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Scn10a CRISPRa kit - CRISPR gene activation of mouse sodium channel, voltage-gated, type X, alpha
Symbol:	Scn10a
UniProt ID:	Q6QIY3
Synonyms:	Nav1.8; PN3; SNS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)