

## Product datasheet for **GA203782**

### Mouse *Scn3a* activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Scn3a CRISPRa kit - CRISPR gene activation of mouse sodium channel, voltage-gated, type III, alpha
Symbol:	Scn3a
UniProt ID:	<a href="#">A2ASI5</a>
Synonyms:	Gm1000; Nav1.3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)