

Product datasheet for **GA203786**

Mouse *Scn9a* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Scn9a CRISPRa kit - CRISPR gene activation of mouse sodium channel, voltage-gated, type IX, alpha
Symbol:	Scn9a
Components:	GA203786G1, Scn9a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA203786G2, Scn9a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA203786G3, Scn9a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q62205
Synonyms:	mKIAA4197; Nav1.7; PN1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)